


---


**JULIET**<sup>®</sup>  
Ti PO/OL

---

BY SPINEART



POSTERIOR  
LUMBAR  
Ti CAGES





# CONTENT

CONCEPT AND DESIGN	<b>PAGE 04</b>
IMPLANTS	<b>PAGE 05</b>
TECHNICAL FEATURES	<b>PAGE 07</b>
INSTRUMENT SET	<b>PAGE 09</b>
INSTRUMENTS	<b>PAGE 10</b>
SURGICAL TECHNIQUE	<b>PAGE 12</b>



## CONCEPT AND DESIGN

In 2006, to accompany the ROMEO® posterior fusion system, Spineart developed a range of interbody devices to achieve 360° fusion: the JULIET interbody system.

Named after William Shakespeare's characters Romeo and Juliet, the two systems complement each other perfectly.

The JULIET<sub>PO</sub>, JULIET<sub>OL</sub>, JULIET<sub>AN</sub> and JULIET<sub>TL</sub> are designed to be used with the ROMEO 2 system for a reliable, efficient and easy-to-use platform to achieve fusion.

Building on the success and experience acquired with our PEEK range, Spineart developed a new Titanium range, featuring the *Ti-LIFE*<sub>Technology</sub>, a state-of-the-art porous, interconnected structure replicating the trabecular bone geometry.

With each product development, Spineart is relentlessly driven by the same motto: Quality, Innovation and Simplicity.

### AT A GLANCE

**Ti-LIFE TECHNOLOGY**

**OSTEOCONDUCTIVE**

**IMAGING PERFORMANCES**

**COMPLETE RANGE**



### INDICATIONS


The JULIET<sub>Ti</sub> range is indicated for arthrodesis of the lumbar spine at one level or two contiguous levels from L2 to S1 in patients with:

- Degenerative pathology, including symptomatic disc degeneration, recurrent hernia, degenerative spondylolisthesis;
- Isthmic spondylolisthesis;

These patients should be skeletally mature and have had six months of non-operative therapy. Supplemental fixation/stabilization as well as additional bone grafting material are required.




# JULIET<sub>Ti PO</sub>



REFERENCES			
LORDOSIS 6°		W	
H07	L24	10.5	JUT-P6 24 07-S
H08	L24	10.5	JUT-P6 24 08-S
H09	L24	10.5	JUT-P6 24 09-S
H10	L24	10.5	JUT-P6 24 10-S
H11	L24	10.5	JUT-P6 24 11-S
H12	L24	10.5	JUT-P6 24 12-S
H13	L24	10.5	JUT-P6 24 13-S
H14	L24	10.5	JUT-P6 24 14-S



JULIET Ti OL



REFERENCES			
LORDOSIS 6°		W	
H07	L28	10.5	JUT-O6 28 07-S
H08	L28	10.5	JUT-O6 28 08-S
H09	L28	10.5	JUT-O6 28 09-S
H10	L28	10.5	JUT-O6 28 10-S
H11	L28	10.5	JUT-O6 28 11-S
H12	L28	10.5	JUT-O6 28 12-S
H13	L28	10.5	JUT-O6 28 13-S
H14	L28	10.5	JUT-O6 28 14-S

REFERENCES			
LORDOSIS 6°		W	
H07	L32	10.5	JUT-O6 32 07-S
H08	L32	10.5	JUT-O6 32 08-S
H09	L32	10.5	JUT-O6 32 09-S
H10	L32	10.5	JUT-O6 32 10-S
H11	L32	10.5	JUT-O6 32 11-S
H12	L32	10.5	JUT-O6 32 12-S
H13	L32	10.5	JUT-O6 32 13-S
H14	L32	10.5	JUT-O6 32 14-S

OPTIONAL REFERENCES*			
LORDOSIS 6°		W	
H07	L36	10.5	JUT-O6 36 07-S
H08	L36	10.5	JUT-O6 36 08-S
H09	L36	10.5	JUT-O6 36 09-S
H10	L36	10.5	JUT-O6 36 10-S
H11	L36	10.5	JUT-O6 36 11-S
H12	L36	10.5	JUT-O6 36 12-S
H13	L36	10.5	JUT-O6 36 13-S
H14	L36	10.5	JUT-O6 36 14-S


\* On special request



# TECHNICAL FEATURES



## Ti-LIFE TECHNOLOGY

- The osteoconductive structure mimics the bone trabecular geometry and is designed to promote bone in-growth. This technology is based on a propriety algorithm associated with a unique additive manufacturing process, commonly referred to as 3D printing.




## SMOOTH BULLET NOSE

- The cages feature a smooth bullet self-distracting nose and polished chamfer. This design is aimed for easy insertion, enabling distraction of the intervertebral space while mitigating the risk of damage to the endplates, nerve roots and soft tissue.



## IMAGING PERFORMANCE

- The JULIET<sub>Ti</sub> design features an overall reduced density to optimize the imaging performance.




*Imaging on specimen*



# TECHNICAL FEATURES


## BONE GRAFT

- In addition to the osteoconductive properties of the *Ti-Life*<sub>Technology</sub>, the large windows allow for an extensive bone graft area. Therefore 100 % of the cage surface is dedicated to bone fusion without compromising the mechanical properties of the cage.



## COMPLETE RANGE

- JULIET<sub>Ti PO/OL</sub> cages are available in a wide range of options, to address different patient anatomies, and various surgical approach techniques. For a detailed list of cages please refer to pages 5 and 6 of this guide.



## STREAMLINED AND COMPACT INSTRUMENTATION

- The Combo instrument set provides a complete, modular and compact solution.




# INSTRUMENT SET

## PO/OL COMBO SET



PREPARATION TRAY

OLIF/PLIF TRAY



### UNIVERSAL CONTAINER

CONTAINER - BASE	JUL-BX 10 01-N
------------------	----------------

### PREPARATION TRAY

#	DESCRIPTION	REFERENCE
	UNIVERSAL TRAY	JUL-BX 10 02-N
	UNIVERSAL RACK	JUL-BX 10 05-N
01	T HANDLE	HAN-SI MD TE-N
02	PADDLE DISTRACTOR	JUL-IN 00 05-N JUL-IN 00 06-N JUL-IN 00 07-N
03	DISC SHAVER	JUL-IN 01 07-N JUL-IN 00 08-N JUL-IN 00 09-N JUL-IN 00 10-N JUL-IN 00 11-N JUL-IN 00 12-N JUL-IN 00 13-N JUL-IN 00 14-N
04	MODULAR STRAIGHT HANDLE	HAN-SI SH ST-N
05	SLAP HAMMER	HAN-SS SH 01-N
06	"3 IN 1" COMPACTION BASE	JUT-IN 00 01-N
07	IMPACTION CAP	HAN-SS SH 02-N
08	NERVE ROOT RETRACTOR	DYN-IP 00 05-N
09	COMPACTOR	JUL-IN 14 00-N

### OLIF/PLIF TRAY

#	DESCRIPTION	REFERENCE
	OLIF/PLIF TRAY	JUL-BX 10 04-N
	OLIF/PLIF RACK	JUL-BX 10 11-N
10	IMPLANT HOLDER	DYN-IP 00 01-N
11	IMPLANT HOLDER MICRO	JUL-IO 00 01-N
12	TRIAL IMPLANT JULIET <sub>Ti PO/OL</sub>	JUT-IN 01 07-N JUT-IN 01 08-N JUT-IN 01 09-N JUT-IN 01 10-N JUT-IN 01 11-N JUT-IN 01 12-N JUT-IN 01 13-N JUT-IN 01 14-N
13	CURETTE	JUL-IN 15 00-N

\*For further information on the JULIET<sub>Ti PO/OL</sub> instrumentation, please refer to the corresponding surgical technique.



# INSTRUMENTS

**T HANDLE**

**HAN-SI MD TE-N**



**PADDLE DISTRACTOR**

**JUL-IN 00 XX-N**



**DISC SHAVERS**

**JUL-IN 0X XX-N**




**MODULAR STRAIGHT HANDLE**

**HAN-SI SH ST-N**



**COMPACTOR**

**JUL-IN 14 00-N**



**SLAP HAMMER**

**HAN-SS SH 01-N**



**CURETTE**

**JUL-IN 15 00-N**



**« 3 IN 1 » COMPACTION BASE**

**JUT-IN 00 01-N**



# INSTRUMENTS

**IMPACTION CAP**

**HAN-SS SH 02-N**



**NERVE ROOT RETRACTOR**

**DYN-IP 00 05-N**



**IMPLANT HOLDER PO/OL**

**DYN-IP 00 01-N**



**IMPLANT HOLDER MICRO**

**JUL-IO 00 01-N**



**TRIAL IMPLANT JULIET**  
Ti PO/OL

**JUT-IN 01 XX-N**



# SURGICAL TECHNIQUE : JULIET<sub>Ti OL</sub>

## STEP 1



### DISCECTOMY AND PREPARATION OF THE ENDPLATES

Partially remove the facet joints. Once the approach is done, distract the disc space with the modular paddle distractor, previously assembled with the modular straight handle, or the T handle, for a better rotation.

Proceed to the discectomy.

Prepare and freshen the endplates using the 1mm increment disc shavers. A curette can also be used.

For an optimal protection of the dura, a nerve root retractor is available.

INSTRUMENT	REFERENCE
NERVE ROOT RETRACTOR	DYN-IP 00 05-N
CURETTE	JUL-IN 15 00-N
DISC SHAVER	JUL-IN 0X XX-N
PADDLE DISTRACTOR	JUL-IN 00 XX-N
T HANDLE	HAN-SI MD TE-N
MODULAR STRAIGHT HANDLE	HAN-SI SH ST-N

## STEP 2



### SELECTION OF THE IMPLANT SIZE

To determine the right cage to implant, it is mandatory to use dedicated OL implant trials.

Each implant trial represents the 3 different lengths.

To insert the implant trials, physicians can connect the impaction cap to the modular straight handle to gently hammer on the assembly.

Once satisfied with the selected trial size, proceed to fluoroscopic controls to confirm the correct sizing.

You can use the slap hammer to remove the implant trial.

**Note:** These implant trials can also be used to further rasp the endplates.

INSTRUMENT	REFERENCE
TRIAL IMPLANT JULIET <sub>Ti PO/OL</sub>	JUT-IN 01 XX-N
NERVE ROOT RETRACTOR	DYN-IP 00 05-N
MODULAR STRAIGHT HANDLE	HAN-SI SH ST-N
SLAP HAMMER	HAN-SS SH 01-N
IMPACTION CAP	HAN-SS SH 02-N



# SURGICAL TECHNIQUE : JULIET<sub>Ti OL</sub>

## STEP 3




### CAGE PREPARATION

Select the corresponding cage.  
Connect it with the implant holder.  
Position the cage on the compaction base. Fill it with bone graft or bone substitute.

INSTRUMENT	REFERENCE
IMPLANT HOLDER	DYN-IP 00 01-N
IMPLANT HOLDER MICRO	JUL-IO 00 01-N
« 3 IN 1 » COMPACTION BASE	JUT-IN 00 01-N
COMPACTOR	JUL-IN 14 00-N

## STEP 4



### INSERTION

Insert the cage into the disc space while protecting the dura with the nerve root retractor.  
It is possible to gently hammer on the implant holder handle to ease the insertion of the implant.  
Note that the placement angle of the JULIET<sub>Ti OL</sub> implant is 40 degrees from the median plane.


INSTRUMENT	REFERENCE
NERVE ROOT RETRACTOR	DYN-IP 00 05-N
IMPLANT HOLDER	DYN-IP 00 01-N
IMPLANT HOLDER MICRO	JUL-IO 00 01-N



# SURGICAL TECHNIQUE : JULIET<sub>Ti OL</sub>

## FINAL CONSTRUCT

---



The JULIET<sub>Ti OL</sub> cages should be used with a supplemental posterior fixation system, as described in the ROMEO<sub>2</sub>, ROMEO<sub>2 MIS</sub>, and ROMEO<sub>2 PAD</sub> surgical techniques, or an anterior fixation system.

Compression forceps should be used for final compression of the construct.



# SURGICAL TECHNIQUE : JULIET<sub>Ti PO</sub>

## STEP 1



### DISCECTOMY AND PREPARATION OF THE ENDPLATES


Remove partially the facet joints. Once the approach is done, distract the disc space thanks to modular paddle distractor, previously assembled with the modular straight handle, or the T handle, for a better rotation.

Proceed to the discectomy.

Prepare and freshen the endplates using the 1mm increment disc shavers. A curette can also be used . For an optimal protection of the dura, a nerve root retractor is available.

INSTRUMENT	REFERENCE
NERVE ROOT RETRACTOR	DYN-IP 00 05-N
CURETTE	JUL-IN 15 00-N
DISC SHAVER	JUL-IN 0X XX-N
PADDLE DISTRACTOR	JUL-IN 00 XX-N
T HANDLE	HAN-SI MD TE-N
MODULAR STRAIGHT HANDLE	HAN-SI SH ST-N

## STEP 2



### SELECTION OF THE IMPLANT SIZE

To determine the right cage to implant, it is mandatory to use dedicated PO implant trials.

To insert the implant trials, physicians can connect the impaction cap to the modular straight handle to gently hammer on the assembly.

Once satisfied with the selected trial size, proceed to fluoroscopic controls to confirm the correct sizing.

You can use the slap hammer to remove the implant trial.

**Note:** These implant trials can also be used to further rasp the endplates.

INSTRUMENT	REFERENCE
TRIAL IMPLANT JULIET <sub>Ti PO/OL</sub>	JUT-IN 01 XX-N
NERVE ROOT RETRACTOR	DYN-IP 00 05-N
MODULAR STRAIGHT HANDLE	HAN-SI SH ST-N
SLAP HAMMER	HAN-SS SH 01-N
IMPACTION CAP	HAN-SS SH 02-N



# SURGICAL TECHNIQUE : JULIET<sub>Ti PO</sub>

## STEP 3



### CAGE PREPARATION


Select the corresponding cage.

Position the cage on the compaction base. Fill it with bone graft or bone substitute.

Once the cage is prepared, connect it with the implant holder.

INSTRUMENT	REFERENCE
IMPLANT HOLDER	DYN-IP 00 01-N
IMPLANT HOLDER MICRO	JUL-IO 00 01-N
« 3 IN 1 » COMPACTION BASE	JUT-IN 00 01-N
COMPACTOR	JUL-IN 14 00-N

## STEP 4



### INSERTION

Insert the cage into the disc space by protecting the dura with the nerve root retractor.

It is possible to gently hammer on the implant holder handle to ease the insertion of the implant.

INSTRUMENT	REFERENCE
NERVE ROOT RETRACTOR	DYN-IP 00 05-N
IMPLANT HOLDER	DYN-IP 00 01-N
IMPLANT HOLDER MICRO	JUL-IO 00 01-N




# SURGICAL TECHNIQUE : JULIET<sub>Ti PO</sub>

## STEP 5

### INSERTION OF THE SECOND CAGE

Repeat steps 3 and 4 to ensure the placement of the second cage.




INSTRUMENT	REFERENCE
IMPLANT HOLDER	DYN-IP 00 01-N
IMPLANT HOLDER MICRO	JUL-IO 00 01-N

## FINAL CONSTRUCT

The JULIET<sub>Ti PO</sub> cages should be used with a supplemental posterior fixation system, as described in the ROMEO 2 surgical technique, or an anterior fixation system.

Physicians should use the compression forceps to compress the final construct.







spineart.com

**SPINEART SA**  
Chemin du Pré-Fleuri 3  
1228 Plan-les-Ouates  
Switzerland

**SPINEART USA**  
9130 Irvine Center Drive, Suite 150  
Irvine, CA 92618  
United States of America



0617-V5 ref. JUT-BR OL 31-E

