

# Femur Retrograde Nail System



# Femur Retrograde Nail System


## Implant Ordering Information

### Femur Retrograde Nail

| Cat No.   | Size       | Cat No.   | Size       |
|-----------|------------|-----------|------------|
| 503-49200 | Φ9x200 mm  | 503-41200 | Φ11x200 mm |
| 503-49250 | Φ9x250 mm  | 503-41250 | Φ11x250 mm |
| 503-49300 | Φ9x300 mm  | 503-41300 | Φ11x300 mm |
| 503-49320 | Φ9x320 mm  | 503-41320 | Φ11x320 mm |
| 503-49340 | Φ9x340 mm  | 503-41340 | Φ11x340 mm |
| 503-49360 | Φ9x360 mm  | 503-41360 | Φ11x360 mm |
| 503-49380 | Φ9x380 mm  | 503-41380 | Φ11x380 mm |
| 503-40200 | Φ10x200 mm | 503-42200 | Φ12x200 mm |
| 503-40250 | Φ10x250 mm | 503-42250 | Φ12x250 mm |
| 503-40300 | Φ10x300 mm | 503-42300 | Φ12x300 mm |
| 503-40320 | Φ10x320 mm | 503-42320 | Φ12x320 mm |
| 503-40340 | Φ10x340 mm | 503-42340 | Φ12x340 mm |
| 503-40360 | Φ10x360 mm | 503-42360 | Φ12x360 mm |
| 503-40380 | Φ10x380 mm | 503-42380 | Φ12x380 mm |

### Nail Locking Screw Ø6.5

| Cat No.   | Size | Cat No.   | Size |
|-----------|------|-----------|------|
| 526-65030 | 30mm | 526-65060 | 60mm |
| 526-65032 | 32mm | 526-65062 | 62mm |
| 526-65034 | 34mm | 526-65064 | 64mm |
| 526-65036 | 36mm | 526-65066 | 66mm |
| 526-65038 | 38mm | 526-65068 | 68mm |
| 526-65040 | 40mm | 526-65070 | 70mm |
| 526-65042 | 42mm | 526-65075 | 75mm |
| 526-65044 | 44mm | 526-65080 | 80mm |
| 526-65046 | 46mm | 526-65085 | 85mm |
| 526-65048 | 48mm | 526-65090 | 90mm |
| 526-65050 | 50mm |           |      |
| 526-65052 | 52mm |           |      |
| 526-65054 | 54mm |           |      |
| 526-65056 | 56mm |           |      |
| 526-65058 | 58mm |           |      |



### Nail Locking Screw Ø5.0

| Cat No.   | Size |
|-----------|------|
| 520-25026 | 26mm |
| 520-25028 | 28mm |
| 520-25030 | 30mm |
| 520-25032 | 32mm |
| 520-25034 | 34mm |
| 520-25036 | 36mm |
| 520-25038 | 38mm |
| 520-25040 | 40mm |
| 520-25042 | 42mm |
| 520-25044 | 44mm |
| 520-25046 | 46mm |
| 520-25048 | 48mm |
| 520-25050 | 50mm |



### End Cap

| Cat No.   | Size |
|-----------|------|
| 530-01000 | 0mm  |
| 530-01005 | 5mm  |
| 530-01010 | 10mm |



# Femur Retrograde Nail System

## Instrument Ordering Information

| Cat No.   | Description                                    | Size         |
|-----------|------------------------------------------------|--------------|
| 910-01001 | Protection Sleeve                              | Φ16.0/13.1mm |
| 910-01002 | Drill Sleeve for Guide Pin                     | Φ3.1mm       |
| 910-01003 | Troca                                          | Φ3.0mm       |
| 910-01004 | Open Reamer                                    | Φ 13.0mm     |
| 910-01005 | Curved Awl                                     | Φ 13.0mm     |
| 910-03038 | Primary Awl                                    |              |
| 910-03051 | Guide Pin                                      | Φ 3.0mm      |
| 910-01212 | Assemble Handle for Retrograde                 | Left         |
| 910-01312 | Assemble Handle for Retrograde                 | Right        |
| 910-01213 | Target Guide for Retrograde                    |              |
| 910-01313 | Target Guide for Retrograde, A/P               |              |
| 910-01142 | Assemble Bolt for Retrograde                   |              |
| 910-01119 | Impact Handle for Retrograde                   |              |
| 901-01040 | Drill Bit                                      | Φ4.0x180mm   |
| 901-01143 | Drill Bit                                      | Φ4.3x280mm   |
| 901-01050 | Drill Bit                                      | Φ5.0x280mm   |
| 910-01043 | Protection Sleeve for Femur/Retro, Holding     |              |
| 910-01244 | Drill Sleeve for Retrograde                    |              |
| 910-01145 | Troca                                          | Φ 5.0mm      |
| 910-01243 | Protection Sleeve for Femur/Retro, Non-Holding |              |
| 910-01020 | Extract Bolt                                   |              |
| 910-01120 | Extract Bolt with Tip                          |              |
| 910-03030 | Assemble/Extract Rod                           |              |
| 910-03060 | Wrench                                         |              |
| 910-03029 | Slotted Mallet                                 |              |
| 910-01021 | Assemble Bolt Driver                           | L Ball       |
| 910-01038 | Assemble Bolt Driver                           | T Ball       |
| 910-01048 | End Cap Driver                                 |              |
| 910-01049 | Coupling Screw For End Cap Driver              |              |
| 910-03226 | Screw Driver Shaft, Holding(Small Chuck)       | 3.5mm        |
| 910-03126 | Screw Driver, Holding                          | 3.5mm        |
| 910-01150 | Protection Sleeve, Short (free hand)           |              |
| 910-01151 | Drill Sleeve, Short (free hand)                | Φ 3.5mm      |
| 910-01152 | Troca, Short (free hand)                       | Φ 3.5mm      |
| 910-01050 | Protection Sleeve, Long (free hand)            |              |
| 910-01051 | Drill Sleeve, Long (free hand)                 | Φ 3.5mm      |
| 910-01052 | Troca, Long (free hand)                        | Φ 3.5mm      |
| 910-01026 | Universal Chuck With T- Handle                 |              |
| 910-03125 | Depth Gauge for Nail Locking Screw             |              |
| 910-90007 | Container for Retrograde Nail                  |              |



Sales Office :

#104, #204, Jungang innotceh B/D , 148, Sagimakgol-ro,  
Joongwon-gu, Seongnam-si, Korea  
Tel. 82-31-732-0631 Fax. 82-31-732-0632

Manufacturer :

#F105-1, #101~104, Sangsan-dong, Gwangju Techno-park,  
958-3, Daechon-dong, Buk-gu, Gwangju-si, Korea  
Tel. 82-62-602-7460 Fax. 82-62-602-7461

[www.tradimedics.com](http://www.tradimedics.com)

