

---

ROMEO® 2

---


BY SPINEART

POLY

AXIAL

PEDICLE

SCREW



DEGENERATIVE TL

# ROMEO<sup>®</sup>2



STERILE  
IMPLANT



## STREAMLINED TIP

The screw tip is designed to allow an effortless and self-centering insertion of the screw.

## POLYAXIAL HEAD

The 50° polyaxiality of the screw head is designed to tolerate alignment discrepancies. Monoaxial and spondylolisthesis screws are also available.

## LOW PROFILE IMPLANTS

The low profile of the ROMEO<sup>®</sup>2 implant head is designed to enable an atraumatic implantation and minimize anatomical interference.

## COMPACT SET

Only one box of specific and intuitive instruments is needed for degenerative cases.

S W I S S M A D E



---

ROMEO® 2

---

BY SPINEART

# POLY AXIAL SCREW

